



# Recombinant Mouse DNA polymerase alpha subunit B (Pola2)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP018301MO                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Uniprot No.</b>       | P33611                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Immunogen Species</b> | Mus musculus (Mouse)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Sequence</b>          | MAVSTQQLAE ELQIFGLDYE DLSLEKLAEL CVLYRQTEDG MVSELI AFCT<br>SAGKTCLTVD ILNSFEYEVL NKKLSKAWHS ASKDSGHAGT RDIVSIQELI<br>EAEETETLL SSYTTPSKGP LKRVSSTPET PLTKRSVAAR SPRQLLSPSS<br>FSPSATPSQK YTSRTNRGEV VTTFGSAQGL SWSGRGGSGS VSLKVVGDPE<br>PLTGSYKAMF QQLMGVREVL TSKIEELGSE LKEHHKIEAF TPLLVPAQEP<br>VILLGQIGCD SNGKLNSKSV ILEGDQEHSY GAQIPVDLSE LKEYSLFPGQ<br>VVIMEGFNTT GRRLTATKLY EGVPLPFYQP TEEEGASEQT MVVVACGPYT<br>TSDSITYDPL LDLIAINRD QPDVCILFGP FLDKHEQVE NCKLTSPFED<br>VFKQCLRTVI EGTRSSGSHL VFVPSLRDVH HEPVYPQPPF TFSELSREDK<br>KRVQFVSEPC SLSINGVMFG LTSTDLLFHI GAEEIFSSSG TSDRFSRVLK<br>HILTQRSYYP LYPPHEDMAI DYENFYTYAQ LPVTPDVFIV PSELYFVKD<br>IFGCVCVNPG RLTKGQVGGT FGRLYLRRQP KAMDGGGRQG LSVAAQVVRI |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Target Names</b>      | Pola2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Protein Names</b>     | Recommended name: DNA polymerase alpha subunit B Alternative name(s): DNA polymerase alpha 70 kDa subunit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Expression Region</b> | 1-600                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Protein Length</b>    | Full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                                                                                                                                                       |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                                                                                                                                                |