



Recombinant Bovine NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 9, mitochondrial (NDUFA9)

Product Code	CSB-EP015634BO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34943
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	LHHAV IPHGKGGRRSS VSGIVATVFG ATGFLGRYVV NHLGRMGSSQV IVPHRCEPYD TMHLRPMGDL GQIIFMDWNG RDKDSIRRAV EHSSVVINLV GREWETQNFDFEDVFKIPQ AIAQVSKEAG VEKFIHISHL NADIKSSSKY LRSAVGEKE VRETFPEATI IKPAEIFGRE DRFLNYFANI RWFGGVPLIS LGKKTVKQPV YIVDVTKGII NAIKDPDARG KTFAFVGPSSR YLLFDLVQYV FAVAHRPFLP YPLPHFAYRW IGRLEISPF EPWTTTRDKVE RIHTTDKILP HPLGLEDLGV EATPLELKAI EVLRRHRTYR WLSSEIEDVQ PAKTIPTSGP
Source	E.coli
Target Names	NDUFA9
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 9, mitochondrial Alternative name(s): Complex I-39kD Short name= CI-39kD NADH-ubiquinone oxidoreductase 39 kDa subunit
Expression Region	36-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.