



Recombinant Schizosaccharomyces pombe Mitochondrial protein import protein mas5 (mas5)

Product Code	CSB-EP524885SXV
Abbreviation	Recombinant SchizosaccharoMyces pombe mas5 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O74752
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MVKETKLYEVLNVDVTASQAELKKAYRKLALKYHPDKNPAGDKFKEISRAYEI LADEEKRATYDRFGEEGLQGGGADGGMSADDLFASFFGGGMFGGGMPRGP RKGKDLVHTIKVTLEDLYRGKTTKLALQKKVICPKCSGRGGKEGSVKSCASCN GSGVKFITRAMGPMIQRMQMTCPCNGAGETIRDEDRCCKECDGAKVISQRKIL TVHVEKGMHNGQKIVFKEEGEQAPGIIPGDVIFVIDQKEHPRFKRSGDHLFYEA HVDLLTALAGGQIVVEHLDDRWLTIIPGECIRPNELKVLPGQGMLSQRHHQP GNLYIRFHVDFPEPNFATPEQLALLEKALPPRKIESAPKNAHTEECVLATVDPT EKVRIDNNVDPTTATSMDEDEDEEGGHPGVQC
Research Area	Others
Source	E.coli
Target Names	mas5
Protein Names	Recommended name: Mitochondrial protein import protein mas5
Expression Region	1-404aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	51.5 kDa
Protein Length	Full Length of Mature Protein
Image	

