



Recombinant Human Perilipin-1 (PLIN1), Escherichia Coli

Product Code	CSB-EP526447HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O60240
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAVNKGLTLL DGDLPSEQENV LQRVLQLPVV SGTCECFQKT YTSTKEAHPL VASVCNAYEK GVQSASSLAA WSMEPVVRRLL STQFTAANEL ACRGLDHLLE KIPALQYPPE KIASSELKDTI STRLRSARNS ISVPIASTSD KVLGAALAGC ELAWGVARDT AEFAANTRAG RLASGGADLA LGSIEKVVEY LLPPDKEESA PAPGHQQAQK SPKAKPSLLS RVGALTNTLS RYTVQTMARA LEQGHTVAMW IPGVVPLSSL AQWGASVAMQ AVSRRRSEVR VPWLHSLAAA QEEDHEDQTD TEGEDTEEEE ELETEENKFS EVAALPGPRG LLGGVAHTLQ KTLQTTISAV TWAPAAVLGM AGRVLHLTPA PAVSSTKGRA MSLSDALKGV TDNVVDTVVH YVPLPRLSLM EPESEFRDID NPPAEVERRE AERRASGAPS AGPEPAPRLA QPRRSLRSAQ SPGAPPGPGL EDEVATPAAP RPGFPAVPRE KPKRRVSDSF FRPSVMEPIL GRTHYSQLRK KS
Source	E.coli
Target Names	PLIN1
Protein Names	Recommended name: Perilipin-1 Alternative name(s): Lipid droplet-associated protein
Expression Region	1-522
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein coats lipid storage droplets in adipocytes, thereby protecting them until they can be broken down by hormone-sensitive lipase. The encoded protein is the major cAMP-dependent protein kinase substrate in adipocytes and, when unphosphorylated, may play a role in the inhibition of lipolysis. Alternatively spliced transcript variants varying in the 5' UTR, but encoding the same protein, have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.