




Human Prealbumin protein

Product Code	CSB-NP003401h
Relevance	Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain. Detected in serum and cerebrospinal fluid (at protein level). Highly expressed in choroid plexus epithelial cells. Detected in retina pigment epithelium and liver.
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Tested Applications	ELISA, WB, SDS-PAGE
Form	Liquid
Storage Buffer	PBS, pH 7.4
Alias	Transthyretin, ATTR, TBPA
Product Type	Native Protein
Sensitivity	Not test
Purity	≥95% (SDS-PAGE)
Sequence	Full length protein
Lead Time	2 weeks
Research Area	Cardiovascular/Stem Cells Lineage/Cancer
Source	Purified from Human serum
Protein Names	Human Prealbumin protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.