




Rabbit Transferrin

Product Code	CSB-NP003801Rb
Relevance	Transferrins are iron binding transport proteins which can bind two Fe ³⁺ ions in association with the binding of an anion, usually bicarbonate. It is responsible for the transport of iron from sites of absorption and heme degradation to those of storage and utilization. Serum transferrin may also have a further role in stimulating cell proliferation.
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Tested Applications	ELISA, WB, SDS-PAGE
Form	Liquid
Storage Buffer	PBS, pH 7.4
Product Type	Native Proteins
Sensitivity	Not test
Purity	≥95% (SDS-PAGE)
Sequence	Full length protein
Lead Time	2 weeks
Research Area	Cardiovascular
Source	Purified from Rabbit Serum
Protein Names	Rabbit Transferrin

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.