




Chicken Ovomuroid

Product Code	CSB-NP004001C
Relevance	N/A
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Tested Applications	ELISA, WB, SDS-PAGE
Form	Liquid
Storage Buffer	PBS, pH 7.4
Product Type	Native Proteins
Sensitivity	Not test
Purity	>95% (SDS-PAGE)
Sequence	Full length protein
Lead Time	2 weeks
Research Area	Cardiovascular
Source	Purified from Chicken egg white
Protein Names	Chicken Ovomuroid

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.