




Ovalbumin

Product Code	CSB-NP004301C
Relevance	Storage protein of egg white. Lack protease inhibitory activity.
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Tested Applications	ELISA, WB, SDS-PAGE
Form	Liquid
Storage Buffer	PBS, pH 7.4
Alias	Allergen Gal d II, Egg albumin, Plakalbumin
Product Type	Native Protein
Sensitivity	Not test
Purity	≥95% (SDS-PAGE)
Sequence	Full length protein
Lead Time	2 weeks
Research Area	Cardiovascular/Stem Cells Lineage/Cancer
Source	Purified from Chicken egg white
Protein Names	Ovalbumin

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.