




Guinea pig Hemoglobin protein

Product Code	CSB-NP004801Gp
Relevance	Involved in oxygen transport from the lung to the various peripheral tissues.
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Tested Applications	ELISA, WB, SDS-PAGE
Form	Liquid
Storage Buffer	PBS, pH 7.4
Alias	N/A
Product Type	Native Protein
Sensitivity	Not test
Purity	≥95% (SDS-PAGE)
Sequence	Full length protein
Lead Time	1 week
Research Area	Cardiovascular
Source	Purified from Guinea pig erythrocytes
Protein Names	Guinea pig Hemoglobin protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.