




Mouse Immunoglobulin G2b

Product Code	CSB-NP008301m
Relevance	N/A
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Tested Applications	ELISA, WB, SDS-PAGE
Form	Liquid
Storage Buffer	PBS, pH 7.4
Alias	N/A
Product Type	Native Protein
Sensitivity	Not test
Purity	≥95% (SDS-PAGE)
Sequence	Full length protein
Research Area	Immunology
Source	Purified from pooled normal Mouse serum
Protein Names	Mouse Immunoglobulin G2b

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.