




Recombinant Human Alpha-L-iduronidase (IDUA)

Product Code	CSB-YP011000HU
Abbreviation	Recombinant Human IDUA protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P35475
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	APHLVHVDAARALWPLRRFWRSTGFCPPLPHSQADQYVLSWDQQLNLAYVG AVPHRGIKQVRTHWLLELVTTTRGSTGRGLSYNFTHLDGYLDLLRENQLLP LMGSASGHFTDFEDKQQVFEWKDLVSSLARRYIGRYGLAHVSKWNFETWNE PDHHDNDVSMQMGFLNYYDACSEGLRAASPALRLGGPGDSFHTPPRSPL SWGLLRHCHDGTNFFTGEAGVRLDYISLHRKRGARSSISILEQEKVVAQQIRQLF PKFADTPIYNDEADPLVGWSLPQPWRADVYAAMVVKVIAQHQNLLLANTTSA FPYALLSNDNAFLSYHPHPFAQRTLRTARFQVNNTRPPHVQLLRKPVLTAMGLL ALLDEEQLWAEVSQAGTVLDSNHTVGVLASAHRPQGPADAWRAAVLIYASDD TRAHPNRSVAVTLRLRGVPPGGLVYVTRYLDNGLCSPDGEWRRRLGRPVFPT AEQFRMRAAEDPVAAAPRPLPAGGRLTLR PALRLPSLLLHVHVCARPEKPPG QVTRLRALPLTQGQLVLVWSDEHVGSKCLWTYEIQFSQDGKAYTPVSRKPST FNLVVFSPDTGAVSGSYRVRALDYWARPGPFSDPVYPYLEVPVPRGPPSPGNP
Research Area	Metabolism
Source	Yeast
Target Names	IDUA
Protein Names	Recommended name: Alpha-L-iduronidase EC= 3.2.1.76
Expression Region	28-653aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	71.9kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.