



Recombinant Mouse EGF-containing fibulin-like extracellular matrix protein 1 (Efemp1)

Product Code	CSB-BP007450MO
Abbreviation	Efemp1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8BPB5
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	QYT EETITYTQCT DGYEWDPIRQ QCKDIDECDI VPDACKGGMK CVNHYGGYLC LPKTAQIIVN NEHPQQETPA AEASSGATTG TVAARSMATS GVVPGGGFMA SATAVAGPEV QTGRNNFVIR RNPADPQRIP SNPSHRIQCA AGYEQSEHNV CQDIDECTSG THNCRDQVC INLRGSFTCQ CLPGYQKRGE QCVDIDECTV PPYCHQRCVN TPGSFYCQCS PGFQLAANNY TCVDINECDA SNQCAQQCYN ILGSFICQCN QGYELSSDRL NCEDIDECRT SSYLCQYQCV NEPGKFSCMC PQGYEVVRSR TCQDINECET TNECREDEMC WNYHGGFRCY PRNPCQDHYV LTSENRCVCP VSNTMCRELP QSIVYKYMSI RSDRSVPSDI FQIQATMIYA NTINTFRIKS GNENGEFYLR QTSPVSAMLV LVKSLSGPRE YIVDLEMLTV SSIGTFRTSS VLRLTIIVGP FSF
Source	Baculovirus
Target Names	Efemp1
Protein Names	Recommended name: EGF-containing fibulin-like extracellular matrix protein 1 Alternative name(s): Fibulin-3 Short name= FIBL-3
Expression Region	18-493
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a member of the fibulin family of extracellular matrix glycoproteins. Like all members of this family, the encoded protein contains tandemly repeated epidermal growth factor-like repeats followed by a C-terminus fibulin-type domain. This gene is upregulated in malignant gliomas and may play a role in the aggressive nature of these tumors. Mutations in this gene are associated with Doyme honeycomb retinal dystrophy. Alternatively spliced transcript variants that encode the same protein have been described.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.