



Recombinant Mouse von Willebrand factor (Vwf), partial

Product Code	CSB-BP025960MO
Abbreviation	Vwf
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8CIZ8
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	TEKPRDRP STARCSLFGD DFINTFDETM YSFAGGCSYL LAGDCQKRSF SILGNFQDGK RMSLSVYLGE FFDIHLFANG TVMQGDQSIG MPYASQGLYL ELEAGYYKLS SETFGFAARI DGNGNFQVLM SDRHFNKTCG LCGDFNIFAE DDFRTQEGTL TSDPYDFANS WALSSSEQRC KRASPPSRNC ESSSGDMHQA MWEQCQLLKT ASVFARCHPL VDPESFVALC EKILCTCATG PECACPVILLE YARTCAQEGM VLYGWTDHSA CRPACPAGME YKECVSPCPR TCQSLSINEV CQQQCVDGCS CPEGELLED RCVQSSDCPC VHAGKRYPPG TSLSQDCNTC ICRNSLWICS NEECPGECLV TGQSHFKSFD NRYFTFSGIC QYLLARDCED HTFSIVIETM QCADDPDAVC TRSVSVRLSA LHNSLVKCLKH GGAVGIDGQD VQLPFLQGD LRIQHTVMASV RLSYAEDLQM DWDGRGRLLV KLSPVYSGKT CGLCGNYNGN KGDDFLTPAG LVEPLVVDFG NAWKLQGDSC DLRRQHSDPC SLNPRLTRFA EEACALLTSS KFEACHHAVS PLPYLQNCRY DVCSCSDSRD CLCNAVANYA AECARKGVHI GWREPGFCAL GCPQGQVYLQ CGNSCNLTCR SLSLPDEECS EVCLEGCYCP PGLYQDERGD CVPKAQCPCY YDGEFQPAD IFSDHHTMCY CEDGFMHCTT SGTLGSLLPD TVLSSPLSHR SKR
Source	Baculovirus
Target Names	Vwf
Protein Names	Recommended name: von Willebrand factor Short name= vWF Cleaved into the following chain: 1. von Willebrand antigen 2 Alternative name(s): von Willebrand antigen II
Expression Region	23-763
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Target Details	The glycoprotein encoded by this gene functions as both an antihemophilic



factor carrier and a platelet-vessel wall mediator in the blood coagulation system. It is crucial to the hemostasis process. Mutations in this gene or deficiencies in this protein result in von Willebrand s disease. An unprocessed pseudogene has been found on chromosome 22.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.