




NDRG2 Antibody

Product Code	CSB-PA891788ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UN36
Immunogen	Recombinant Human Protein NDRG2 protein (152-371AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein NDRG2 (N-myc downstream-regulated gene 2 protein) (Protein Syld709613) [Cleaved into: Protein NDRG2, N-terminally processed], NDRG2, KIAA1248 SYLD
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NDRG2

Image



Western blot

All lanes: NDRG2 antibody at 4.93µg/ml + A431 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 41, 40, 37, 38 kDa

Observed band size: 41 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.