




PSMB10 Antibody

Product Code	CSB-PA018877ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40306
Immunogen	Recombinant Human Proteasome subunit beta type-10 protein (1-273AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome subunit beta type-10 (EC 3.4.25.1) (Low molecular mass protein 10) (Macropain subunit MECL-1) (Multicatalytic endopeptidase complex subunit MECL-1) (Proteasome MECL-1) (Proteasome subunit beta-2i), PSMB10, LMP10 MECL1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PSMB10

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA018877ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.