





PEX19 Antibody

Product Code	CSB-PA017802DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40855
Immunogen	Recombinant Human Peroxisomal biogenesis factor 19 protein (1-299AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisomal biogenesis factor 19 (33 kDa housekeeping protein) (Peroxin-19) (Peroxisomal farnesylated protein), PEX19, HK33 PXF
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PEX19

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA017802DSR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA017802DSR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.