





TFCP2 Antibody

Product Code	CSB-PA621517DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12800
Immunogen	Recombinant Human Alpha-globin transcription factor CP2 protein (1-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-globin transcription factor CP2 (SAA3 enhancer factor) (Transcription factor LSF), TFCP2, LSF SEF
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TFCP2

Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA621517DSR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human thymus tissue using CSB-PA621517DSR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.