





GSTZ1 Antibody

Product Code	CSB-PA009994DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43708
Immunogen	Recombinant Human Maleylacetoacetate isomerase protein (1-216AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Maleylacetoacetate isomerase (MAAI) (EC 5.2.1.2) (GSTZ1-1) (Glutathione S-transferase zeta 1) (EC 2.5.1.18), GSTZ1, MAAI
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GSTZ1

Image



Western blot
 All lanes: GSTZ1 antibody at 0.41µg/ml + Mouse liver tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 25, 18, 20 kDa
 Observed band size: 25 kDa



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA009994DSR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.