





CDC7 Antibody

Product Code	CSB-PA005023ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00311
Immunogen	Recombinant Human Cell division cycle 7-related protein kinase protein (275-574AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell division cycle 7-related protein kinase (CDC7-related kinase) (HsCdc7) (huCdc7) (EC 2.7.11.1), CDC7, CDC7L1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CDC7

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA005023ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA005023ESR2HU at dilution of 1:100