



CRYGD Antibody

Product Code	CSB-PA006020DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07320
Immunogen	Recombinant Human Gamma-crystallin D protein (1-174AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-crystallin D (Gamma-D-crystallin) (Gamma-crystallin 4), CRYGD, CRYG4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CRYGD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.