





SERPIND1 Antibody

Product Code	CSB-PA021080ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05546
Immunogen	Recombinant Human Heparin cofactor 2 protein (20-290AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heparin cofactor 2 (Heparin cofactor II) (HC-II) (Protease inhibitor leuserpin-2) (HLS2) (Serpine D1), SERPIND1, HCF2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SERPIND1

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA021080ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human plasma tissue using CSB-PA021080ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.