





LAMTOR5 Antibody

Product Code	CSB-PA010161ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43504
Immunogen	Recombinant Human Ragulator complex protein LAMTOR5 protein (1-91AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ragulator complex protein LAMTOR5 (Hepatitis B virus X-interacting protein) (HBV X-interacting protein) (HBX-interacting protein) (Late endosomal/lysosomal adaptor and MAPK and MTOR activator 5), LAMTOR5, HBXIP XIP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LAMTOR5

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA010161ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA010161ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.