




IL31RA Antibody

Product Code	CSB-PA854162ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NI17
Immunogen	Recombinant Human Interleukin-31 receptor subunit alpha protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-31 receptor subunit alpha (IL-31 receptor subunit alpha) (IL-31R subunit alpha) (IL-31R-alpha) (IL-31RA) (Cytokine receptor-like 3) (GLM-R) (hGLM-R) (Gp130-like monocyte receptor) (Gp130-like receptor) (ZcytoR17), IL31RA, CRL3 GPL
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL31RA

Image



Western blot
 All lanes: IL31RA antibody at 1.45µg/ml + Mouse kidney tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 83, 87, 75, 42, 78, 71, 59, 67, 64, 70, 85 kDa
 Observed band size: 83 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.