




TCF19 Antibody

Product Code	CSB-PA023295ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y242
Immunogen	Recombinant Human Transcription factor 19 protein (1-290AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor 19 (TCF-19) (Transcription factor SC1), TCF19, SC1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TCF19

Image



Western blot

All lanes: TCF19 antibody at 0.96µg/ml

Lane 1: 293T whole cell lysate


Lane 2: Rat thymus tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 37 kDa

Observed band size: 37 kDa



Immunohistochemistry of paraffin-embedded

human testis tissue using CSB-

PA023295ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.