





TCIM Antibody

Product Code	CSB-PA882079ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NR00
Immunogen	Recombinant Human Transcriptional and immune response regulator protein (1-106AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcriptional and immune response regulator (Thyroid cancer protein 1) (TC-1), C8orf4 TC1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TCIM

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA882079ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human lung cancer tissue using CSB-PA882079ESR1HU at dilution of 1:100