





FAIM2 Antibody

Product Code	CSB-PA866298DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWQ8
Immunogen	Recombinant Human Protein lifeguard 2 protein (1-105AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein lifeguard 2 (Fas apoptotic inhibitory molecule 2) (Neural membrane protein 35) (Transmembrane BAX inhibitor motif-containing protein 2), FAIM2, KIAA0950 LFG LFG2 NMP35 TMBIM2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAIM2

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA866298DSR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA866298DSR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.