





LY86 Antibody

Product Code	CSB-PA013252ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95711
Immunogen	Recombinant Human Lymphocyte antigen 86 protein (21-162AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lymphocyte antigen 86 (Ly-86) (Protein MD-1), LY86, MD1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LY86

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA013252ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA013252ESR1HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.