




NOL3 Antibody

Product Code	CSB-PA015921DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60936
Immunogen	Recombinant Human Nucleolar protein 3 protein (1-208AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleolar protein 3 (Apoptosis repressor with CARD) (Muscle-enriched cytoplasmic protein) (Myp) (Nucleolar protein of 30 kDa) (Nop30), NOL3, ARC NOP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NOL3

Image


Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA015921DSR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.