




YEATS4 Antibody

Product Code	CSB-PA026252ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95619
Immunogen	Recombinant Human YEATS domain-containing protein 4 protein (1-227AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	YEATS domain-containing protein 4 (Glioma-amplified sequence 41) (Gas41) (NuMA-binding protein 1) (NuBI-1) (NuBI1), YEATS4, GAS41
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	YEATS4

Image



Western blot

All lanes: YEATS4 antibody at 1.44µg/ml

Lane 1: A549 whole cell lysate

Lane 2: Rat gonad tissue

Lane 3: U251 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 26 kDa

Observed band size: 26 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.