




ZC3HC1 Antibody

Product Code	CSB-PA801246ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86WB0
Immunogen	Recombinant Human Nuclear-interacting partner of ALK protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear-interacting partner of ALK (Nuclear-interacting partner of anaplastic lymphoma kinase) (hNIPA) (Zinc finger C3HC-type protein 1), ZC3HC1, NIPA
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ZC3HC1

Image



Western blot

All lanes: ZC3HC1 antibody at 1.08µg/ml + HeLa whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 56, 54, 48 kDa

Observed band size: 54 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.