




# DNAJA2 Antibody

<b>Product Code</b>	CSB-PA007000DSR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O60884
<b>Immunogen</b>	Recombinant Human DnaJ homolog subfamily A member 2 protein (233-412AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	DnaJ homolog subfamily A member 2 (Cell cycle progression restoration gene 3 protein) (Dnj3) (Dj3) (HIRA-interacting protein 4) (Renal carcinoma antigen NY-REN-14), DNAJA2, CPR3 HIRIP4
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	DNAJA2

## Image



### Western blot

All lanes: DNAJA2 antibody at 3.79µg/ml

Lane 1: A375 whole cell lysate

Lane 2: PC-3 whole cell lysate

Lane 3: HepG2 whole cell lysate

Lane 4: Hela whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 46 kDa

Observed band size: 46 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.