





INPP5J Antibody

Product Code	CSB-PA620891ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15735
Immunogen	Recombinant Human Phosphatidylinositol 4,5-bisphosphate 5-phosphatase A protein (60-360AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositol 4,5-bisphosphate 5-phosphatase A (EC 3.1.3.56) (Inositol polyphosphate 5-phosphatase J), INPP5J, PIB5PA PIPP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	INPP5J

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA620891ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA620891ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.