





MYLK3 Antibody

Product Code	CSB-PA658362ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q32MK0
Immunogen	Recombinant Human Myosin light chain kinase 3 protein (1-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin light chain kinase 3 (EC 2.7.11.18) (Cardiac-MyBP-C-associated Ca/CaM kinase) (Cardiac-MLCK), MYLK3, MLCK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYLK3

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA658362ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA658362ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.