



# ANKRD36B Antibody

<b>Product Code</b>	CSB-PA818689ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8N2N9
<b>Immunogen</b>	Recombinant Human Ankyrin repeat domain-containing protein 36B protein (1-230AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Ankyrin repeat domain-containing protein 36B (CLL-associated antigen KW-1), ANKRD36B
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Others
<b>Target Names</b>	ANKRD36B
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.