




NRBF2 Antibody

Product Code	CSB-PA822209DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96F24
Immunogen	Recombinant Human Nuclear receptor-binding factor 2 protein (10-220AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear receptor-binding factor 2 (NRBF-2) (Comodulator of PPAR and RXR), NRBF2, COPR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NRBF2

Image



Western blot

All lanes: NRBF2 antibody at 5.51µg/ml

Lane 1: HepG2 whole cell lysate


Lane 2: A549 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 33, 27, 32 kDa

Observed band size: 33 kDa



Immunohistochemistry of paraffin-embedded

human placenta tissue using CSB-

PA822209DSR2HU at dilution of 1:100

