




SKAP2 Antibody

Product Code	CSB-PA021355DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75563
Immunogen	Recombinant Human Src kinase-associated phosphoprotein 2 protein (1-120AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Src kinase-associated phosphoprotein 2 (Pyk2/RAFTK-associated protein) (Retinoic acid-induced protein 70) (SKAP55 homolog) (SKAP-55HOM) (SKAP-HOM) (Src family-associated phosphoprotein 2) (Src kinase-associated phosphoprotein 55-related protein) (Src-associated adapter protein with PH and SH3 domains), SKAP2, PRAP RA70 SAPS SCAP2 SKAP55R
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SKAP2

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA021355DSR2HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.