





TIMM17A Antibody

Product Code	CSB-PA860770ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99595
Immunogen	Recombinant Human Mitochondrial import inner membrane translocase subunit Tim17-A protein (1-171AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial import inner membrane translocase subunit Tim17-A (Inner membrane preprotein translocase Tim17a), TIMM17A, MIMT17 TIM17 TIM17A TIMM17
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TIMM17A

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA860770ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA860770ESR1HU at dilution of 1:100

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.