




TPP2 Antibody

Product Code	CSB-PA024114ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29144
Immunogen	Recombinant Human Tripeptidyl-peptidase 2 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripeptidyl-peptidase 2 (TPP-2) (EC 3.4.14.10) (Tripeptidyl aminopeptidase) (Tripeptidyl-peptidase II) (TPP-II), TPP2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TPP2

Image



Western blot

All lanes: TPP2antibody at 3.42µg/ml

Lane 1: Hela whole cell lysate

Lane 2: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 138 kDa

Observed band size: 138 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA024114ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.