





CHIA Antibody

Product Code	CSB-PA874858ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZP6
Immunogen	Recombinant Human Acidic mammalian chitinase protein (1-368AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acidic mammalian chitinase (AMCase) (EC 3.2.1.14) (Lung-specific protein TSA1902), CHIA
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CHIA


Image



Western blot
 All lanes: CHIA antibody at 0.74µg/ml + Mouse stomach tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 53, 41, 34 kDa
 Observed band size: 53 kDa



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA874858ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA874858ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.