



FGF14 Antibody

Product Code	CSB-PA008620ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92915
Immunogen	Recombinant Human Fibroblast growth factor 14 protein (1-252AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor 14 (FGF-14) (Fibroblast growth factor homologous factor 4) (FHF-4), FGF14, FHF4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FGF14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.