





FRZB Antibody

Product Code	CSB-PA835696ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92765
Immunogen	Recombinant Human Secreted frizzled-related protein 3 protein (166-325AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Secreted frizzled-related protein 3 (sFRP-3) (Frezzed) (Fritz) (Frizzled-related protein 1) (FrzB-1), FRZB, FIZ FRE FRP FRZB1 SFRP3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FRZB

Image



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA835696ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA835696ESR1HU at dilution of 1:100