




HERPUD1 Antibody

Product Code	CSB-PA622984ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15011
Immunogen	Recombinant Human Homocysteine-responsive endoplasmic reticulum-resident ubiquitin-like domain member 1 protein (1-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homocysteine-responsive endoplasmic reticulum-resident ubiquitin-like domain member 1 protein (Methyl methanesulfonate (MMF)-inducible fragment protein 1), HERPUD1, HERP KIAA0025 MIF1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	HERPUD1

Image



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA622984ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.