





CNMD Antibody

Product Code	CSB-PA012854ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75829
Immunogen	Recombinant Human Leukocyte cell-derived chemotaxin 1 protein (75-334AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leukocyte cell-derived chemotaxin 1 (Chondromodulin) [Cleaved into: Chondrosurfactant protein (CH-SP); Chondromodulin-1 (Chondromodulin-I) (ChM-I)], CNMD, CHMI LECT1 MYETS1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CNMD

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA012854ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA012854ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.