





LEFTY2 Antibody

Product Code	CSB-PA012858ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00292
Immunogen	Recombinant Human Left-right determination factor 2 protein (57-366AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Left-right determination factor 2 (Endometrial bleeding-associated factor) (Left-right determination factor A) (Protein lefty-2) (Protein lefty-A) (Transforming growth factor beta-4) (TGF-beta-4), LEFTY2, EBAF LEFTA LEFTYA TGFB4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LEFTY2


Image



Western blot
 All lanes: LEFTY2 antibody at 2µg/ml + A431 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 41, 38 kDa
 Observed band size: 41 kDa



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA012858ESR2HU at dilution of 1: 100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA012858ESR2HU at dilution of 1: 100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.