




TAGLN2 Antibody

Product Code	CSB-PA023106DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37802
Immunogen	Recombinant Human Transgelin-2 protein (1-199AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transgelin-2 (Epididymis tissue protein Li 7e) (SM22-alpha homolog), TAGLN2, KIAA0120
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	TAGLN2

Image



Western blot

All lanes: TAGLN2 antibody at 4.04µg/ml

Lane 1: 293T whole cell lysate

Lane 2: K562 whole cell lysate


Lane 3: MCF-7 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 23, 25 kDa

Observed band size: 23 kDa



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA023106DSR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.