





ZDHHC17 Antibody

Product Code	CSB-PA811600ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IUH5
Immunogen	Recombinant Human Palmitoyltransferase ZDHHC17 protein (170-310AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Palmitoyltransferase ZDHHC17 (EC 2.3.1.225) (Huntingtin yeast partner H) (Huntingtin-interacting protein 14) (HIP-14) (Huntingtin-interacting protein 3) (HIP-3) (Huntingtin-interacting protein H) (Putative MAPK-activating protein PM11) (Putative NF-kappa-B-activating protein 205) (Zinc finger DHHC domain-containing protein 17) (DHHC-17), ZDHHC17, HIP14 HIP3 HYPH KIAA0946
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ZDHHC17

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA811600ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA811600ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.