





CDKN2D Antibody

Product Code	CSB-PA005095DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55273
Immunogen	Recombinant Human Cyclin-dependent kinase 4 inhibitor D protein (1-166AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cyclin-dependent kinase 4 inhibitor D (p19-INK4d), CDKN2D
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CDKN2D

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA005095DSR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA005095DSR1HU at dilution of 1:100