




# COL9A3 Antibody

<b>Product Code</b>	CSB-PA005759DSR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q14050
<b>Immunogen</b>	Recombinant Human Collagen alpha-3(IX) chain protein (445-684AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Collagen alpha-3(IX) chain, COL9A3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	COL9A3

## Image



### Western blot

All lanes: COL9A3antibody at 1.36μg/ml

Lane 1: MCF-7 whole cell lysate


Lane 2: A549 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 64 kDa

Observed band size: 64 kDa



### Immunohistochemistry of paraffin-embedded

human skin tissue using CSB-

PA005759DSR1HU at dilution of 1:100