




DGKA Antibody

Product Code	CSB-PA006832ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23743
Immunogen	Recombinant Human Diacylglycerol kinase alpha protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diacylglycerol kinase alpha (DAG kinase alpha) (EC 2.7.1.107) (80 kDa diacylglycerol kinase) (Diglyceride kinase alpha) (DGK-alpha), DGKA, DAGK DAGK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DGKA

Image



Western blot

All lanes: DGKA antibody at 6.48µg/ml

Lane 1: Hela whole cell lysate

Lane 2: Jurkat whole cell lysate


Lane 3: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 83, 14, 29 kDa

Observed band size: 83 kDa



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA006832ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.