




FGF19 Antibody

Product Code	CSB-PA008624ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95750
Immunogen	Synthesized peptide derived from human FGF-19
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:3000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FGF19

Image



Western Blot

Positive WB detected in: 293T whole cell lysate, HeLa whole cell lysate, HepG2 whole cell lysate, MCF7 whole cell lysate

All lanes: FGF19 Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 25 kDa

Observed band size: 25 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.